



Touchscreen-Display:

Easy readout and adjustment of the process parameters



Front-panel with LCD-Touchscreen



Transducer-System

Hand held unit or for stand-clamping



Optional Accessories



Pedal-switch



Stainless steel stand
for clamping the transducer-
system

Sonotrodes

Type	Length	Front surface diameter (mm)	Material	Volumen ml/max
MR1 – 30 Tip 1	approx. 80	1	Spec. Titanium	0.1 – 0.5
MR2 – 30 Tip 2	approx. 80	2	Spec. Titanium	2 - 50
MR3 – 30 Tip 3	approx. 80	3	Spec. Titanium	5 - 100
MR7 – 30 Tip 7	approx. 80	7	Spec. Titanium	10 - 250
MR7 – 30 Tip 7 With gasket	approx. 80	7	Spec. Titanium	10 - 250
MR7 – 30L Tip 7 With gasket	approx. 160	7	Spec. Titanium	10 - 250
MR10 – 30 Tip 10 (only LS – MR100)	approx. 80	10	Spec. Titanium	20 - 500